

Designation: E 1259 - 9401

Standard Test Method for Evaluation of Antimicrobials in <u>Distillate Liquid</u> Fuels (Based on Preliminary Screening and Compatibility) Boiling Below 390°C¹

This standard is issued under the fixed designation E 1259; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ϵ) indicates an editorial change since the last revision or reapproval.

1. Scope

- 1.1 This test method is designed to evaluate antimicrobial agents for the prevention of microbial-<u>ly</u> indfluenced deterioration of <u>distillate liquid</u> fuels (as defined by Specification <u>D 396 as fuel</u>) or <u>D 396, D 910, D 975, D 1655, D 2069, D 2880, D 3699, D 4818 and D 6227), system deterioration, or both. Note 1—A knowledge</u>
 - 1.2 Knowledge of microbiological techniques is required for these procedures.
- 1.23 It is the responsibility of the investigator to determine whether Good Laboratory Practice (GLP) is required and to follow them where appropriate (40 CFR, 160), or as revised.
- 1.34 This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use. See caution statement, Note 2, in Section 8.

2. Referenced Documents

2.1 ASTM Standards:

¹ This test method is under the jurisdiction of ASTM Committee E-35 on Pesticides and is the direct responsibility of Subcommittee E35.15 on Antimicrobial Agents. Current edition approved June 15, 1994. April 10, 2001. Published August 1994. June 2001. Originally published as E 1259 – 88. Last previous edition E 1259 – 8894.



D 396 Specification for Fuel Oils²

D-4054 Practice 910 Specification for Evaluating the Compatibility of Additives with Aviation-Turbine Aviation Gasolines²

D 975 Specification for Diesel Fuel Oils²

D 1655 Specification for Aviation Turbine Fuels²

D 2069 Specification for Marine Fuels²

D 2880 Specification for Gas Turbine Fuels³

D 3699 Specification for Kerosine³

D 4814 Specification for Automotive Spark-Ignition Engine Fuel⁴

D 6227 Specification for Grade 82 Unleaded Aviation Gasoline⁴

<u>D 6469 Guide to Microbial Contamination in Fuels</u> and <u>Aircraft</u> Fuel <u>System Materials</u> <u>Systems</u>⁴

2.2 Federal Standards:

² Annual Book of ASTM Standards, Vol 05.01.



- 40 CFR, Part 160, Good Laboratory Practice Standards⁴ 79, Fuels and Fuel Additives Registration Regulations⁵
- 40 CFR, Part 152, Pesticide Registration and Classification Procedures⁵

3. Terminology

- 3.1 Definitions of Terms Specific to This Standard:
- 3.1.1 antimicrobials, n—chemical agents that are cidal or static to microorganisms.
- 3.1.2 *microbial-induced deterioration*, *n*—decomposition/degradation of material (fuel) or making unsuitable for use, as a result of metabolic activity or the presence of microbes
- 3.1.3 *microbicide*, *n*—an agent that kills microbes: bacterial vegetative cells, fungal vegetative cells and spores, algae and protozoa. This term is applied to chemical agents that kill microbes.

4. Summary of Test Method

34.1 This test method is conducted on a reference fuel—for, and determininges the antimicrobial efficacy under well-defined conditions that include specific inocula *Pseudomonas aeruginosa*, American Type Culture Collection, (ATCC) No. 33988, *Hormoconis resinae*, ATCC No. 20495, and *Yarrowia tropicalis* (formerly *Candida tropicalis*, ATCC No. 18138; as well as water/fuel ratios, and time of containment. It is designed for destructive sampling at regular intervals during bottom water buildup. This test method allows for impact of fuel/water partitioning and time, on the antimicrobial agent, as well as the effect of continual rechallenge. Every 2 weeks, water phase is increased by 0.25 % while concomitantly, a paired system is destructively tested. Thus, at 4 weeks, there is an equivalent increase in the water phase to 0.5 %, at 6 weeks 0.75 %, and at 8 weeks 1.0 %. At each sampling time interval, treated and untreated aliquots are checked for the three types of organisms in the initial inoculum. These counts are coupled with gross observations of each system for biofilm formation and interfacial growth.

4. Significance and Use

4.1 The procedure should be used to evaluate the relative efficacy of microbicides in distillate fuels. The effect of environmental conditions including a variety of fuel additives, metal surfaces, and climatology are variables that can be included in specific tests using this protocol.

5. Significance and Use

- 5.1 Guide D 6469 details the types of problems associated with uncontrolled microbial growth in fuels and fuel systems. Treatment with effective antimicrobial agents is one element of contamination control strategy.
- 5.2 The procedure should be used to evaluate the relative efficacy of microbicides in distillate fuels. The effect of environmental conditions, such as a variety of fuel additives, metal surfaces, and climatology, are variables that can be included in specific tests using this protocol.
- 5.3 This method addresses product performance issues only. Regulatory Agencies restrict and control the use of both pesticides (in the U.S.:40 CFR 152) and fuel additives (40 CFR 79). Regardless of performance in this method, antimicrobials must only be used in compliance with applicable regulations. Specific industries, for example, the aviation industry, may place further restrictions on chemicals used for fuel treatment.

6. Apparatus

- 56.1 Colony Counter—Any of several types, for example, a Quebec Colony Counter may be used.
- 56.2 Incubator—Any incubator capable of maintaining temperature of 30 to $35 \pm 2^{\circ}$ C may be used.
- 56.3 Sterilizer—Any suitable steam sterilizer capable of producing the conditions of sterility is acceptable.
- 56.4 Separatory Funnels—Eight—Sixteen 1-L funnels.
- 56.5 Ring Stand, suitable for supporting separatory funnel.
- 56.6 Vortex—Mixer.

67. Reagents and Materials

- 67.1 Petri Dishes—100 by 15 mm required for performing standard plate count.
- 67.2 Bacteriological Pipets—10.0 mL and 1.1, or 2.2 mL capacity.
- 67.3 Water Dilution Bottles—Any sterilizable glass container having a 150–200 mL capacity and tight closure may be used. 6.4 Distillate Fuel.⁵

³ Annual Book of ASTM Standards, Vol 05.02.

⁴ Available from the Superintendent

⁴ Annual Book of Documents, U.S. Government Printing Office, Washington, DC 20402. ASTM Standards, Vol 05.03.

⁵ I-H CAT diesel fuel is available

⁵ Available from Howell Hydrocarbons Inc., San Antonio, TX 78223-3531: the Superintendent of Documents, U.S. Government Printing Office, Washington, DC 20402.

6.5 Synthetic Bottom Water.

7.4 Fuel.⁶

67.5 Synthetic Bottom Water.⁷

7.6 Soy Peptone Casein Digest Agar.

67.7 Sabouraud Dextrose Agar.

67.8 Agar, Bacteriological Grade.

67.9 Potassium Tellurite Solution—sterile 1 %.

67.10 Gentamicin Sulfate—50 µg/mL.

67.11 Plate Count Agar.

67.12 Potato Dextrose Agar.

78. Inoculum

78.1 <u>Inoculum Preparation and Maintenance:</u>

8.1.1 *Inoculum Revitalization*—Cultures are *Pseudomonas aeruginosa*, ATCC No. 33988, *Hormoconis resinae*, ATCC No. 20495, and *Yarrowia tropicalis* (formerly *Candida tropicalis*), ATCC No. 18138. Obtain lyophilized preparations from ATCC. Before initiating fuel antimicrobial tests, revitalize each of the three cultures in accordance with the instructions contained with each culture.

8.1.2 Maintenance and Preparation of Inocula—All three cultures are transferred from-specific agar slants of a specified agar, (a) Pseudomonas aeruginosa, plate count agar; (Plate Count Agar), (b) Hormoconis resinae Potato Dextrose Agar), and (c) Yarrowia tropicalis, p (Potato-d Dextrose-a Agar) to synthetic bottom water medium in a suitable size screw-cap glass bottle (French square), and then overlaid with 10 times the volume of fuel. This two-phase system is kept at room temperature (18-24°C) for seven days, and the interface with half the bottom water is transferred weekly to a similar system weekly until used. The bacterial levels expected are about 10⁷ CFU/mL, Colony Forming Units, CFU/mL, the yeast levels 10⁶ CFU/mL, and mold levels 10⁴ spores/mL. For the test inoculum, the bacteria are diluted 1:100 while the latter two yeast and molds are diluted 1:10. The counting of the inoculum-(zero time) is done directly from the prepared synthetic bottom water mixture at time zero, just prior to adding inoculum to each-setup at zero time setup, and at each subsequent time point.

Note 1—Caution: In the distillate fuel industry, additives, including biocides, are calculated on a weight per weight basis so that the specific gravity of both the fuel and the biocide (if a liquid formulation) must be taken into account.

89. Procedure

89.1 *The Setup*—For each biocide, four setups are needed plus four for the control. A typical evaluation could include three levels of the biocide in addition to the control system with no biocide.

$$4 \times 3l + 4c = 16ts \tag{1}$$

where:

l = level,

= control, and

ts = total setups.

89.1.1 To each-1_1-L-t_funnel, add 800 mL of test fuel-with and the appropriate level of-biocide.

Note 2—Caution: In the distillate fuel industry, additives, including biocides, are calculated on a weight per weight basis stock biocide so that the specific gravity desired concentration of both the fuel and the biocide-(if is achieved.

9.1.2 Stock biocide is made by adding it to a liquid formulation) must be taken into account.

8.1.2 Next, solvent that will dissolve the biocide.

9.1.3 Next, add 0.25 % 2 mL of fuel volume (2 mL) of synthetic bottom water medium (0.25 % of fuel volume) containing the appropriate levels of each organism. Close the funnel. Shake vigorously for 10 s. Place in ring stand. Open stopper slightly to allow volatile gas to escape. Close and leave for 2 two weeks.

89.2 Sampling—At each time point after 0 weeks, at 2, 4, 6, and 8 weeks, the following protocol is observed. The bottom water fraction including the fuel/water interface and a minimal amount of fuel is bled from the separatory funnel and mixed vigorously in a vortex for 10 s. Before settling, aliquots are removed for plate counts for the yeast, bacteria, yeast and molds with dilutions into mold plate counts. Dilutions, 1:100, 1:1000, and 1:10000, are prepared in sterile synthetic bottom water solution as follows: 1.0 mL into 99 mL for 1:100 and from this dilution subsequent dilutions for 1:1 000 and 1:10 000 are made. water. These will be used for pour platings for bacteria and yeast, respectively. In addition, three by two 0.1 mL portions will be used for bacteria and yeast pour plates, and for spread plates for mold counts; f. For estimating Pseudomonas aeruginosa, use soy casein digest agar;

⁶ Items 6.5-6.12

⁶ Representative fuel samples from each product grade are available from-a variety of media manufacturers and chemical supply companies. all petroleum refiners.

Items 7.5-7.12 are available from a variety of media manufacturers and chemical supply companies.



- for Yarrowia tropicalis, use Sabouraud's-d Dextrose-a Agar with gentamicin 0.5 μg/mL; and for Hormoconis resinae, use 0.01 % potassium tellurite in 1.5 % bacteriological agar.
- 9.3 Reinoculation—Immediately after the sacrificing (destructive sampling) of one-system biocide level, 2 mL of synthetic bottom water solution containing additional inoculum is added to each of the three remaining setups/biocide level. Again, the procedure described in 8.1.2 and 9.1.3 is repeated. After two more weeks, 9.2 is again repeated and the process is repeated at weeks 6 and 8.

10. Results

- 10.1 Comparison of Test and Control—At each interval, the microbiological counts for treated systems will be compared with those of the untreated systems. In addition, gross observation of the condition of each system will be made with the intent of using these data as part of the evaluation.
- Note—3—_2—Yarrowia readily outgrows Hormoconis in—s_Sabouraud making distinction of both groups difficult, if not impossible. Hormoconis resinae is able to grow in a simple, unsupplemented agar, albeit slowly (about 5 days incubation with a tellurite reduction as an indicator of growth). Under these minimal nutritional conditions, the potassium tellurite may also be inhibitory to the yeast.
- 8.3 Reinoculation—Immediately after the sacrificing (destructive sampling) of one-system biocide level, 2 mL of synthetic bottom water solution containing additional inoculum is added to each of the three remaining setups/biocide level. Again, the procedure of shaking, etc. described in 8.1 is repeated. After 2 more weeks, 8.2 is repeated and the process described here is repeated at Weeks 6 and 8.

9. Results

9.1 Comparison of Test and Control—At each interval, the microbiological counts for treated systems will be compared with those of the untreated systems. In addition, gross observation of the condition of each system will be made with the intent of using these data as part of the evaluation.

10. Precision and Bias

10.1 It is not practical to specify the precision of the procedure in Test Method E 1259 because detection and enumeration of microorganisms is subjective and not absolute. Since there is no accepted reference material suitable for the procedure in Test Method E 1259, bias has not been determined.

11. Precision and Bias

11.1 It is not practical to specify the precision of the procedure in Test Method E 1259 because detection and enumeration of microorganisms is subjective and not absolute. Since there is no accepted reference material suitable for the procedure in Test Method E 1259, bias has not been determined.

12. Keywords

1+2.1 antimicrobials; aviation fuels; biodeterioration; diesel; distillate fuels; gasoline; gas-turbine fuels; marine fuels; microbially-induced deterioration

APPENDIX

(Nonmandatory Information)

X1. ALTERNATIVE PROCEDURES

- X1.1 Setups using 1-L French squares or other suitable container, can substitute for separatory funnels. Although the interface is not so readily discerned, the bottom water can be removed by careful pipetting.
- X1.2 Fuel—Potentially, the inclusion of additives may not only increase the growth rate of the inoculum but may also affect the biocide negatively as well as positively. efficacy of the biocide. In any specific testing for fuel suppliers, it would be unrealistic to exclude these additives from the study (see also Practice D 4054).
- X1.3 Bottom Water Level and Time of Storage—The test described here does not consider either long term storage with minimal bottom water (for example, 8 weeks at 0.25 % water) or shock dosing a heavily contaminated system with a water soluble biocide, or bleed off of bottom water with loss of biocide.
- X1.4 *Corrosion*—There is provision for use of metal coupons for both evaluation of biocide corrosivity and microbial-induced corrosion.



X1.5 *Inoculum*—For specific testing, it may be advisable to use contaminated fuel as an inoculum after determination of the identity of the inoculum. contaminant.

ASTM International takes no position respecting the validity of any patent rights asserted in connection with any item mentioned in this standard. Users of this standard are expressly advised that determination of the validity of any such patent rights, and the risk of infringement of such rights, are entirely their own responsibility.

This standard is subject to revision at any time by the responsible technical committee and must be reviewed every five years and if not revised, either reapproved or withdrawn. Your comments are invited either for revision of this standard or for additional standards and should be addressed to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend. If you feel that your comments have not received a fair hearing you should make your views known to the ASTM Committee on Standards, at the address shown below.

This standard is copyrighted by ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, PA 19428-2959, United States. Individual reprints (single or multiple copies) of this standard may be obtained by contacting ASTM at the above address or at 610-832-9585 (phone), 610-832-9555 (fax), or service@astm.org (e-mail); or through the ASTM website (www.astm.org).